

AM2 SYSTEM

MEDICAL SYSTEMS

An Outstanding Value in Orthopedic Radiography



AM2 System Provides Value

- Ideal for private office or clinic, this economical system provides all of the benefits of the AM1 system except the elevating table feature.
- This system also includes a handsome four-way float-top table with a precision tabletop bearing system. Cross-table, weight-bearing knees, chest procedures, table work -- all are possible with the full-featured rotational tubestand.
- Also included is the accurately counterbalanced orthopedic wallstand.
- As with the elevator table system, the AM2 four-way system has electric locks on all table and tubestand motions.



Complement either Americomp orthopedic system with your choice from a complete line of high-quality, state-of-the-art generators. See your Americomp dealer for accompanying literature.