

MEDICAL SYSTEMS

The Americomp AM4 Radiographic System



- This affordable system features a sturdy table that supports loads up to 300 pounds. An extra-long pedestal base adds stability and full bucky travel.
- With a detached floor-to-wall (or to ceiling) tubestand, the AM4 System is capable of a wide range of radiographic procedures including table-top extremities, skull, and abdominal table work.
- When used with the Americomp Wallstand, the system capabilities include standing limbs and chest exams.
- This versatile film receptor may include a stationary fine-line grid or a reciprocating bucky. The cabinet accepts a film cassette tray which accommodates up to a 14" x 17" cassette in either direction.



This photo demonstrates the versatility of the Americomp tubestand.