



Mobile X-ray Imaging is Really on the Move.

Now, X-Factor Thinking lets you take the DRX detector to the field with this rugged, all-in-one solution. It includes the wireless detector, electronics, wireless access point, and a tablet PC—all securely packed in a sturdy, portable carrying case. The DRX-Transportable System is ideal for applications such as the military, veterinary medicine, or travel between nursing homes and home-care patients. You'll have all the performance of a DR imaging room—wherever you need it.



ADVANTAGES AT A GLANCE

- Compact, protective carrying case for the Field Portable kit makes it easy to transport to remote locations.
- Quickly and cost-effectively convert existing analog mobile systems to DR.
- Easy to use, minimal training necessary.
- Wireless capability improves productivity and workflow.
- Available with either DRX-1 or more dose sensitive DRX-1C detector.

CARESTREAM DRX-Transportable System/Field Portable

The Power of DR Imaging, More Portable Than Ever Before.

- Capture and view images immediately at the point of patient care
- Lightweight wireless DRX detector is easy to use in confined spaces
- Self-powered and rechargeable
- Water and dust resistant rugged carrying case for easy transport



X-Factor Thinking

The X-Factor makes the most of your capital spending with configurable, modular components powered by the portable DRX detector. Use the detector on your retrofitted mobile units on your morning rounds... and later, take it to a busy DRX-Evolution room for peak patient volumes. When the X-Factor lets you scale your solution to fit your workflow and budget, the possibilities are endless.

Trusted Partner – CARESTREAM Smart Services

At Carestream, we'll be there when you need us. Our medical imaging solutions are backed by our worldwide support team of product and industry specialists—experts who work closely with you to keep your equipment up and running at peak performance. By delivering a comprehensive scope of service offerings, we provide complete peace of mind over the life of your equipment, so you can focus on what's most important—your patients.

Product Specifications

Key Components:

- Thin electronics box 454 x 368 x 48 mm (17.9 x 14.5 x 1.9 in.) containing the wireless access point, ethernet switch, generator interface board, and power supply
- Electronics box uses the same type of battery as the DRX detector
 - Battery charges while plugged in
- Motion Computing™ tablet PC
 - I7 Processor / 12.1 in. touch-screen / 1.6 kg (3.6 lbs)
 - Rugged design that meets MIL-STD-810G and IP52
 - Dual Hot Swappable battery
- Mounting hardware to mount tablet PC to system
- Carrying case 546 x 457 x 241 mm (21.5 x 18.0 x 9.5 in.)



www.carestream.com

