

3546/3547



.....3546/3547 Vertical Wall Holder.....



3546/3547 Vertical Wall Holder

Designed specifically for upright radiography, the 3546/3547 Vertical Wall Holder meets the requirements of hospitals, clinics, ambulatory care centers, HMO's, group practice and private offices. Extensive vertical bucky travel enables the user with rapid positioning of the patient. The bucky lowers to floor level to permit weight bearing knees and ankle exams.

Combine the 3546/3547 Vertical Wall Holder with one of our many system packages, for the ultimate high performance radiographic system.

System Features

- ▲ Floor-to-wall mounted vertical counterbalanced for increased esthetics and stability.
- ▲ Vertical mechanical lock maintains cabinet height.
- ▲ Fully counterbalanced wall cabinet positions easily to desired height.
- ▲ Dual counterweight cables for added patient/user safety.
- ▲ 17" x 17" L/F 800 reciprocating bucky with radiolucent front panel.
- ▲ Supplied with cassette size sensing tray, and PBL interface.
- ▲ Accepts 14" x 17" and all smaller size cassettes to adapt to any technique.
- ▲ 103 line, 8:1 ratio fine line grid provided for greater film quality.
- ▲ 61.3" of vertical travel permits a wide variety of radiographic exams.
- ▲ AEC compatible.

Comprehensive Equipment Warranties

Specifications subject to change without notice.
 Conforms to CDRH Standards at time of shipment.
 Baked on Durable Medical White Finish.
 Available from stock for quick, simple installations.

