

The Piloter Series is currently awaiting FDA Clearance and cannot be sold or demonstrated in the USA.

Piloter Series

*A Smart Tablet Diagnostic Ultrasound System
Affordable Excellence for All*



With a 13.3 inch LED Monitor and Touch Screen, the Piloter is a compact Color Doppler Ultrasound System for all Diagnostic Applications.

Piloter is a specialized ultrasound product, dedicated for all diagnostic ultrasound exams. With its extremely lightweight case, smart design and built-in customized workflow, Piloter provides a more efficient and convenient experience for General Imaging, MSK, Pain Management, Vascular Imaging, Vascular Access, OB/GYN and much more. This product was built on the HOLO Platform, has 64-beam processing and achieves up to 5000 frames per second. Piloter also has available the wiNeedle optional technology which provides Intelligent Needle Visualization. Some other optional features are the wiLearn educational software and the the transducer expansion adapter. Piloter has a genius design that can be used anywhere.

wisonic
Your needs we care


The Piloter Series is currently awaiting FDA Clearance and cannot be sold or demonstrated in the USA.

Piloter Series

A Smart Tablet Ultrasound System

Affordable Excellence for All

Designed for the demands of
diagnostic ultrasound imaging



The Piloter Series with its fully digital beamformer, supports most applications, including Radiology, Vascular, OB/GYN, Breast Imaging, Pediatrics, MSK, Cardiology, Emergency Medicine, Vascular Access, Small Parts, Shared Service, Superficial Imaging and much more. Some of the features the Piloter Series supports are 2D, Colorflow, PW Doppler, CW Capability for Cardiac, wiNeedle for automatic needle visualization during biopsies on linear transducers, Dual Live Colorflow imaging, Auto Doppler measurements, extended field of view on linear and convex transducers, quick Boot up times, custom user presets and much more. The Piloter has both Ethernet and WIFI connections and fully supports DICOM 3.0.

Some of the Piloter Series Transducers



C5-2B Probe



L15-4NB Linear



LH15-6
Hockystick Linear



EV10-4
Endocavity



P4-1
Phased Cardiac



L10-5 Linear



C8-3
Micro-Convex



Probe Adapter

For information contact **WisonicUSA** at: Sales@WisonicUSA.net or 804-690-3475 or 407-973-4541