

RT100

Four-way Float Top Table

This sturdy table is ideally suited for a variety of general radiographic needs. Its compact size makes the RT100 easy to use in small hospitals, clinics, ambulatory care centers and private offices. The four-way floating table top glides freely on a precision roller bearing system and can support patients up to 700 lbs (318 kg)

Key Highlights

Four-way Float Top

- Glides effortlessly in all directions

Heavy-duty Electric Table Locks with Foot Treadle Release

700 lbs (318 kg) Patient Weight Capacity

All Steel Structure Welded Base

- Provides maximum strength and reliability

Pedestal Base Can House Generator Electronics to Save Space

Smooth, Stain Resistant Fiber Resin Top




Technical Specifications

* All specifications subject to change without notice

RT100 Four-Way Float Top Table

Load capacity:	700 lbs (318 kg)
Table base length:	43 in (109 cm)
Table base width:	29 in (74 cm)
Tabletop movement:	Longitudinal, 51.50 in (131 cm) Transverse, 9.5 in (24 cm)
Tabletop-to-film distance:	3.25 in (8 cm)
Total bucky travel:	Longitudinal 11.25 in (28 cm)
Tabletop length:	86.66 in (220 cm)
Tabletop width:	32 in (81 cm)
Weight:	387 lbs (176 kg)
Power requirements:	120/230-240 VAC
Options:	78 in (198 cm) tabletop length / Variety of grid options / Lateral cassette holder / Rotating cassette tray / Table-side patient handgrips




Rotating cassette tray



Table-side patient handgrips



Lateral cassette holder



DEL MEDICAL



#9000-RT100 01.22 Rev. C

New York

28 Calvert Street
Harrison, NY 10528
(Tel) 800.261.9808
(Fax) 914.835.6111

Illinois

241 Covington Drive
Bloomington, IL 60108
(Tel) 800.800.6006
(Fax) 847.288.7011

www.delmedical.com