



# CFPM400 & CFPM401 Imaging Tables



Integrated Headrest Table

Oakworks Imaging Tables have a large height range helping to create safe patient transfers while allowing doctors to stand comfortably. Our large lifting column is very rigid, ensuring rock solid support for fast imaging and procedures. Available with an integrated headrest or rectangle top.

## Motions

- Motion 1 - Height Range 26" - 44"
- Motion 2 - Trendelenburg Tilt +15°/-12°
- Motion 3 - Longitudinal Travel 10"
- Motion 4 - Lateral Tilt  $\pm 15^\circ$  (CFPM400)
- Motion 4 - Lateral Travel 8" (CFPM401)

## Features

- Choice of integrated headrest or rectangular carbon fiber tops with 56" of metal free imaging
- 2" Radiolucent Comfort Foam<sup>®</sup> table top pad
- Warranty is 3 years parts and 2 years labor

# SPECIFICATIONS

- Large Adjustable Height Range 26" - 44"
- 56" Metal Free Imaging Space
- Cost Effective 4 Movement Table
- 500 lb Patient Weight Capacity
- Warranty: 3 years parts, 2 years labor

**CFPM400** #81772 Integrated #81808 Rectangular **CFPM401** #81802 Integrated #81809 Rectangular

## INCLUDED FEATURES

**Carbon Fiber Top** 56" Metal free imaging space  
**Integrated Headrest** 22" (56cm) wide x 84" (214cm) long  
**Rectangular** 24" (61cm) wide x 80" (203cm) long

## CFPM400 MOVEMENTS

Motion 1	Electric Height Range 26" - 44" (66-112cm) (4" casters w/o pad)
Motion 2	Electric Trendelenburg tilt +15°/-12°
Motion 3	Electric Longitudinal travel 10" (25cm)
Motion 4	Electric Lateral tilt ±15°

## CFPM401 MOVEMENTS

Motion 1	Electric Height Range 26" - 44" (66-112cm) (4" Casters w/o pad)
Motion 2	Electric Trendelenburg tilt +15°/-12°
Motion 3	Electric Longitudinal travel 10" (25cm)
Motion 4	Electric Lateral travel ±4" (10cm)

Controls	All functions hand control
Weight Capacity	500 lb (227kg)
Product Weight	475 lb (216kg)
Included Features	2" (5cm) Radiolucent Comfort Foam™ table top pad, 4" (10cm) locking casters, Crescent face rest pad (included with integrated headrest top), patient restraint straps (2)
Warranty	3 years parts, 2 years labor
Safety Listings	FDA listed, CE marked, cETLus listed, CB Certificate

## OPTIONS *must be ordered with table.*

		#
5" (13cm) casters	Recommended for mobile service providers	65427
Foot Control (CFPM400)	All functions control	75077
Foot Control (CFPM401)	All functions control	75078
Non-USA electrical systems	Please inquire	--

## ACCESSORIES *order any time.*

		#
Arm Hammock	Use with Integrated Headrest Top	48751
Carbon Fiber Arm Board (6" wide)	180° Articulation	83575
Fluoro Extender	Polycarbonate	33975
Interventional Supine Pillow System	Contoured pad and support wedge pad	81821
Spine Positioning System II	8 piece system	81819
T-rail Clamp	6" (15cm) long rail	63336



Rectangular Top Table

# ACCESSORIES

## SPINE POSITIONING SYSTEM II

Radiolucency & adjustability for ideal cervical, thoracic & lumbar visualization. Minimize unwanted movement & reduce procedure time.



## ARM HAMMOCK

This unique product will comfortably & easily support your patient's arms when prone. Easily attached to the Integrated Head Rest table top.



## SINGLE CARBON FIBER ARM BOARD

Fully radiolucent with a 180° range of motion, the arm board requires no mounting or attachment hardware and simply slides under the table top pad. It can be used on any fluoroscopy table & will support arms and feet in a variety of positioning configurations.



## FLUORO EXTENDER

Made of radiolucent, polycarbonate material, it's organic shape is conducive to a variety of table placement applications. 30" long x 23.5" wide x 1/4" thick



## 6" ADJUSTABLE T-RAIL CLAMP

The adjustable T-rail clamp can be positioned along the grooved table top edge. Clamp is also removable.



## 5" CASTER (OPTIONAL)

Recommended for mobile service providers.



PMSB4901

Contact Oakworks  
 717.759.3146  
[www.oakworksmed.com](http://www.oakworksmed.com)  
[medsales@oakworks.com](mailto:medsales@oakworks.com)